

Figure	Microscope, Wavelength	Activation Protocol	Total Activation Time
<i>Figure 1b</i>	LSM880, 488nm	Every 20 sec	1 min 40 sec
<i>Figure 1d</i>	Spinning Disc, 488nm	Global (1000ms) every 20 sec	1 min
<i>Figure 3a-b</i>	LSM880, 488nm	Every 20 sec	4min
<i>Figure 4a-b</i>	LSM880, 488nm	Every 20 sec	4 min
<i>Figure 5</i>	LSM880, 488nm	Every 20 sec	4 min
<i>Figure 6</i>	LSM880, 488nm	Every 20 sec	1 min 20 sec
<i>Figure 7a-b</i>	LSM880, 488nm	Every 20 sec	4 min