

250/membrane top

Cntnap2+/+

Cntnap2-/-

CASPR2@150

cut @100

cut @20

PV @12

10

membrane bottom

cut @100

75

50

GAPDH @37

25

cut @20

100-250

10-20

20-100

Old (errored) GAPDH membrane

250/membrane top

Cntnap2+/+

Cntnap2-/-

CASPR2@150

cut @100

cut @20

PV @12

10

membrane bottom

cut @100

75

50

GAPDH @37

25

cut @20

100-250

10-20

20-100

Old (errored) GAPDH membrane

3 min