



Figure 2 supplement 1. Axotomy in the DNC enables accurate axon growth following the original path

(A) DA9 regeneration 48hr after axotomy in the DNC. The lightning mark indicates the cutting site. The scale bar = 10 μ m.

(B) DA9 regeneration 48hr after axotomy in the dorsal-ventral midline of the commissure. The lightning mark indicates the cutting site. The arrowhead indicates the growth cone. The scale bar = 10 μ m.

(C) Quantification of regenerated DA9 axons following the original axon path or an alternative path. The number indicates number of animals.