



Figure 5—Figure supplement 2. A split-Gal4 driver expressed specifically in MB neurons.

Anti-GFP (green) and anti-Brp (magenta) staining is shown for maximal projections of adult brain and ventral nerve cord expressing *20×UAS-IVS-CD8-GFP* under the control of a mushroom body-specific split-Gal4 driver (MB004B). Inset shows maximal projections of a partial z-stack generated from posterior sections of a different brain, indicating broad expression in MB neuron somata. Scale bars, 100 μm and 50 μm (inset).