




Anti-p190-glued



Input (5 ng)	
GST	
GST-Kaz C	
GST-Kaz C Nt (aa1-176)	
GST-Kaz C Ct (aa250-327)	
GST	
GST-Kaz C	
GST-Kaz C Nt (aa1-176)	
GST-Kaz C Ct (aa250-327)	
Plus dynactin	Minus dynactin



Input (5 ng)	
GST	
GST-Kaz C	
GST-Kaz C Nt (aa1-176)	
GST-Kaz C Ct (aa250-327)	
GST	
GST-Kaz C	
GST-Kaz C Nt (aa1-176)	
GST-Kaz C Ct (aa250-327)	
Plus dynactin	Minus dynactin



Input (5 ng)

GST

GST-Kaz C

GST-Kaz C Nt (aa1-176)

GST-Kaz C Ct (aa250-327)

GST

GST-Kaz C

GST-Kaz C Nt (aa1-176)

GST-Kaz C Ct (aa250-327)

Plus dynactin

Minus dynactin